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
SOP : Standard Operating Procedure

MANUAL OF HISTOPATHOLOGY, CYTOLOGY, AUTOPSY AND MOLECULAR BIOLOGY

2nd Edition

Department of Pathology
Queen Elizabeth Hospital
Hong Kong

Edited by
John K.C. Chan

Contributed by:

John K.C. Chan, S.K. Chan, L.C. Ho, Polly W.Y. Lam, K. Lee,
Ivy S.C. Luk, W.M. Pong, William Y.W. Tsang, C.H. Tse, S.K. Wan

January 1998

44 - Immunohistochemistry Principles

IMMUNOHISTOCHEMISTRY

BASIC PRINCIPLES

Various substances (e.g. protein, hormone) may be specifically localized in tissue (tissue sections or cytologic smears) by immunohistochemical methods. The major techniques whereby the antigen-antibody reaction may be visualized are:

1. Immunofluorescence

- Direct technique: Antibody directly tagged with fluorescent dye (such as fluorescein isothiocyanate, abbreviated as FITC).
- This method is used mostly on frozen sections (such as renal or skin biopsy) because fixed tissues commonly give nonspecific fluorescence.
- Immunofluorescence study, being a relatively simple and short procedure, is performed manually at QEH.

2. Immunoperoxidase or immunoalkaline phosphatase

- Peroxidase or alkaline phosphatase is used to visualize the antigen-antibody reaction, via a bridge detection system. The former oxidizes a substrate (most commonly DAB, which is being used at QEH) into a dark brown colored product. The latter converts the substrate (Fast Red) into a bright red colored product.
- At QEH, Immunoperoxidase and immunoalkaline phosphatase studies are performed using automated immunostainers: either the VENTANA machine (mostly performing immunoperoxidase staining on tissue sections) or the TECHMATE machine (which utilizes the capillary gap method, and mostly used in our laboratory for performing immunoalkaline phosphatase staining on cytologic smears).
- The method currently used in our laboratory is the labeled streptavidin-biotin peroxidase (LSAB-P) or labeled streptavidin-biotin alkaline phosphatase (LSAB-AP) system for the VENTANA automatic immunostainer. The former is the standard method used on tissue sections, performed at 37°C. The latter is the standard method used on cytologic smears (very thinly spread smears or cytospins); it is applied on tissue sections only on special request, such as tissue with abundant endogenous brown pigment.

- The sequence of reagents is as follow for LSAB-P (with buffer wash between steps):

- a) Inhibitor (hydrogen peroxide and sodium azide to inhibit endogenous peroxidase) for 4 minutes.
- b) Primary antibody (mouse origin for monoclonal antibody, and rabbit origin for antiserum) for 32 minutes.
- c) Universal biotinylated secondary antibody (containing a mixture of biotinylated anti-mouse and biotinylated anti-rabbit antibodies in protein stabilizer and preservative) for 8 minutes.
- d) Labeled streptavidin-peroxidase for 8 minutes.
- e) Substrate: DAB-hydrogen peroxide solution for 8 minutes, followed by copper sulfate enhancement for 4 minutes.

- The sequence of reagents is as follow for LSAB-AP (with buffer wash between the steps):

- a) Primary antibody (mouse origin for monoclonal antibody, and rabbit origin for antiserum) for 32

44 - Immunohistochemistry Principles

IMMUNOHISTOCHEMISTRY

BASIC PRINCIPLES

Various substances (e.g. protein, hormone) may be specifically localized in tissue (tissue sections or cytologic smears) by immunohistochemical methods. The major techniques whereby the antigen-antibody reaction may be visualized are:

1. Immunofluorescence

- Direct technique: Antibody directly tagged with fluorescent dye (such as fluorescein isothiocyanate, abbreviated as FITC).
- This method is used mostly on frozen sections (such as renal or skin biopsy) because fixed tissues commonly give nonspecific fluorescence.
- Immunofluorescence study, being a relatively simple and short procedure, is performed manually at QEH.

2. Immunoperoxidase or immunoalkaline phosphatase

- Peroxidase or alkaline phosphatase is used to visualize the antigen-antibody reaction, via a bridge detection system. The former oxidizes a substrate (most commonly DAB, which is being used at QEH) into a dark brown colored product. The latter converts the substrate (Fast Red) into a bright red colored product.
- At QEH, Immunoperoxidase and immunoalkaline phosphatase studies are performed using automated immunostainers: either the VENTANA machine (mostly performing immunoperoxidase staining on tissue sections) or the TECHMATE machine (which utilizes the capillary gap method, and mostly used in our laboratory for performing immunoalkaline phosphatase staining on cytologic smears).
- The method currently used in our laboratory is the labeled streptavidin-biotin peroxidase (LSAB-P) or labeled streptavidin-biotin alkaline phosphatase (LSAB-AP) system for the VENTANA automatic immunostainer. The former is the standard method used on tissue sections, performed at 37°C. The latter is the standard method used on cytologic smears (very thinly spread smears or cytospins); it is applied on tissue sections only on special request, such as tissue with abundant endogenous brown pigment.

- The sequence of reagents is as follow for LSAB-P (with buffer wash between steps):


- a) Inhibitor (hydrogen peroxide and sodium azide to inhibit endogenous peroxidase) for 4 minutes.
- b) Primary antibody (mouse origin for monoclonal antibody, and rabbit origin for antiserum) for 32 minutes.
- c) Universal biotinylated secondary antibody (containing a mixture of biotinylated anti-mouse and biotinylated anti-rabbit antibodies in protein stabilizer and preservative) for 8 minutes.
- d) Labeled streptavidin-peroxidase for 8 minutes.
- e) Substrate: DAB-hydrogen peroxide solution for 8 minutes, followed by copper sulfate enhancement for 4 minutes.

- The sequence of reagents is as follow for LSAB-AP (with buffer wash between the steps):

- a) Primary antibody (mouse origin for monoclonal antibody, and rabbit origin for antiserum) for 32

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




CD3

antigen list		北海道大学医学部附属病院病理部 抗体データベース		検索	全てを 対象
antibody	CD30(Ki-1, Ber-H2)	clone	Ber-H2	isotype	mono
retrieval	pressure cooker	☐ amplification	メーカー	DAKO	
			保管場所	4℃	
濃度	20	QE antibody	DAKO		
data	20-40	control	tonsil		
staining pattern	Only cell membrane or Golgi staining is considered true positive. Diffuse cytoplasmic staining should be disregarded.				
染色性	Directed against the Ki-1 antigen. Normally BerH2 (CD30) stains isolated large lymphoid cells around the lymphoid follicles. Activated lymphoid cells are also CD30 positive, e.g. infectious mononucleosis, Toxoplasmic lymphadenitis, Kikuchi's lymphadenitis. Plasma cells can also be stained. CD30抗原は、活性化B細胞および活性化T細胞に発現しますが、非活性化リンパ球または単球ではみられません。また、CD30はHodgkinリンパ腫のReed-Sternberg細胞に発現します。 Ber-H2 stains the following tumors: (i) Anaplastic large cell lymphoma ("Ki-1 lymphoma") (ii) Some conventional large cell lymphomas (iii) Isolated large cells in some peripheral T-cell lymphoma and B-cell lymphoma (iv) Reed-Sternberg cells of classical Hodgkin's disease, but it does not stain the L&H cells of nodular lymphocyte predominance Hodgkin's disease (v) Lymphomatoid papulosis, and the group of primary cutaneous CD30+ lymphoproliferative disorders (vi) Plasmacytoma (some cases) (vii) Embryonal carcinoma				
備考	ドット状の染まりがgolgi staining. ALCL T/nullにはALKも染めること。				
				価格	39000
				購入日	2000.10.01

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©2005